# Manual

HL888HA

---

**Blodtryksapparat** til overarmen  
**Blood Pressure Monitor** for the upper arm

# Indholdsfortegnelse

---

## **Introduktion**

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---

## **Introduction**

- What is blood pressure? 2
- WHO classification 2
- Variations in blood pressure 3
- Why is it useful to measure blood pressure at home? 4

## **Device description**

- Description of the blood pressure monitor 5
- Description of display symbols 6

## **Correct use**

- Important before use of the unit 8

## **How to use the product**

- Placing the cuff 9
- Changing batteries 10
- Setting Month, Date and Time 10
- How to measure 11
- Storing and recalling date 12
- Maintenance and warning 13

# Introduction

---


## What is blood pressure?

Blood pressure is the pressure or the force of blood flowing against the walls of the arteries. Blood pressure is constantly changing during the course of the cardiac cycle. Systolic blood pressure is the highest pressure in the cycle.

Diastolic blood pressure is the lowest pressure in the cycle. Both pressure readings, the Systolic rate and the Diastolic rate, are information for the physician to evaluate the status of patient's blood pressure. Many factors such as physical activities, anxieties, time of day, can influence one's blood pressure. Blood pressure is typically low in the morning and increases from afternoon to the evening. It is lower in the summer and higher in the winter.

## WHO klassifikation.

Standards for assessing high or low blood pressure without regard to age have been established by the World Health Organization (WHO) as shown in the chart.


---

## Variations in blood pressure

Individual blood pressures vary greatly on a daily and seasonal basis. These variations are even more pronounced in hypertension patients. Normally the blood pressure rises while at work and is at its lowest during sleep. (hypertension: means a person who has high blood pressure symptom.)

The graph below illustrates the variations in blood pressure over a whole day with measurement taken every five minutes. The thick line represents sleep.

The rise in blood pressure at 4 PM ( A in the graph) and 12 PM ( B in the graph) correspond to an attack of pain.

---

### **Why is it useful to measure blood pressure at home?**

Because having one's blood pressure measured by a doctor in a hospital, or clinic tends to stimulate nervousness in a person and may cause the blood pressure to rise. Also, blood pressure changes in different conditions, so judgement on the basis of single measurement is not reliable.

The blood pressure measured first thing in the morning after getting up, before taking any food is known as the fundamental blood pressure.


# Device description

---


Measures human's Systolic blood pressure, Diastolic blood pressure and heart rate on upper arm by oscillometric measurement method.

The product is recommended for use by people over the age of 18 and is for home use and not clinical use.

1. LCD display
2. Air plug
3. Arm cuff
4. Air jack
5. Start/Stop button
6. Recall memory button
7. Set button
8. Mode button
9. 1.5V x 4 alkaline Batteries (AA)

## Description of display symbols

Display mark	Condition/Cause	Corrective action
 <p><b>pulse mark</b></p>	<p>Flashes on the LCD screen when pulse is detected during measurement.</p>	<p>Measurements in progress. Do not move or talk.</p>
 <p><b>replace batteries</b></p>	<p>Appears when the battery voltage is excessively low or the positions of batteries are incorrect.</p>	<p>Replace all four batteries with new ones. Insert the batteries in correct polarities.</p>
 <p><b>measuring error</b></p>	<p>Appears when the accurate blood pressure could not be obtained.</p>	<p>Press "START/STOP" button again to re-measure. Or check cuff if wrapped properly around the arm and according to instructions. Check palm if exerting effort. Check if there is movement during measurement. Check if posture is correct.</p>

# Important before use of the unit





---

- Please rest for at least 5-10 minutes before taking the measurement.
- To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3-5 minutes between measurements. You may need to adjust the wait time according to your personal physiological situation.
- If you have one of the circulatory problems as arteriosclerosis, diabetes, liver disease, kidney disease, severe hypertension, peripheral circulation....., please consult your doctor or healthcare professional before using the devices.
- Wait 30 ~ 45 minutes before measurement if you've just consumed caffeinated beverages or smoked cigarettes.
- Blood pressure measurements should be interpreted by a physician or trained health professional who is familiar with one's medical history. By using the unit regularly and recording the results for the physician to interpret, one can keep physician informed of the continuing trends in one's blood pressure.
- Wrap the cuff snugly around your arm, the cuff must be at the same level as your heart.
- Do not move the unit during measurement, or accurate measurement can not be achieved.
- Take measurement in a relaxed position without talking.
- Do not wrap the cuff over jacket / sweater sleeve , or measurement cannot be done.
- Keep in mind that blood pressure naturally varies from time to time through out the day and is affected by lots of different factors such as smoking, alcohol consumption, medication, and physical activities.
- Blood pressure measurement determined by the unit are equivalent to those obtained by a trained observer using the cuff/ stethoscope auscultation method, and is within the limits prescribed by the Standard of EN 1060-4.

# How to use the product

---

- Please rest for at least 5-10 minutes before taking the measurement.
- Slide the end of arm cuff furthest from the tube through the metal ring to a loop. The smooth cloth should be on the inside of the cuff.
- Put left arm through the cuff loop.
- The bottom of the cuff should be approx. 0.5 inch above elbow.
- Pull the cuff so that the top and bottom edges are tightened around your arm.
- When the cuff is positioned properly, press the Velcro firmly against the pile side of the cuff.
- Sit on a chair and place your arm on the table so that the cuff is at the same level as your heart.
- Relax your arm and turn your arm upward.
- Make sure there are no kinks in the air tube. To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3-5 minutes between measurements. You may need to adjust the wait time according to your personal physiological situation.

## IMPORTANT!


It is extremely important that the cuff is at the same level as the heart and the cuff fits your arm.

## Changing batteries

- Remove the battery cover at the rear side of the unit.
- Insert ALKALINE batteries into the battery compartment, Don't mix new and used batteries
- Make sure the polarities + and - are observed.
- Information will appear on the screen in 3 seconds

## Setting Month, Date and Time

- Press "MODE" button ("month" starts flashing.)
- Press "SET" button to set the correct month (1,2,.....12)
- Press "MODE" button again ("date" starts flashing )
- Press "SET" button to set the correct date
- Press "MODE" button again ("hour" starts flashing )
- Press "SET" button to set the correct hour in 12-hour format.
- Press "MODE" button again ("minute" starts flashing )
- Press "SET" button to set the correct minute (0,1,2,3,....59)
- Press "MODE" button again to switch to normal time. (month:10 , date:5, hour:12 and minute:00 as shown on the picture to the right)

# How to use the product

---

- Wrap the cuff around the arm (refer to "How to use the product")
- Sit upright on the chair to have correct posture. (refer to "How to use the product")
- Press START/STOP button, all elements will display '888' for 3 seconds. Afterwards, the updated time will be shown.
- The monitor will automatically inflate to the level that is right for you.

## **Do not move or talk while taking blood pressure measurement.**

- After the air pressure is increased, it will slowly decrease. When the pulse is detected, the HEART RATE symbol will start flashing.
- The inflation will go to 160 and test if it is enough, if your blood pressure is higher, it will continue.
- After taking the blood pressure measurement, the systolic rate, diastolic rate and pulse rate will be shown on the display for 1 minute, before the unit will automatically switch off.

### **Note!**

**Press START/STOP button anytime to stop measuring.**

---

## Storing and recalling data.

### Storing data:

- After each blood pressure measurement, the Systolic pressure, Diastolic pressure, Pulse rate and the time and date of specific day will be automatically stored.
- The memory holds the latest 99 measurement data. If more than 99 measurements, the memory will automatically clear out the earliest data.

### Recalling data:

- Pressing the "MEMORY" button at normal time, nothing will appear on the display if there is no data in the memory. If there is data in the memory, the latest recorded set of measurement which are the Systolic pressure, Diastolic pressure, Pulse rate, including the recording time and date, will appear for 1 minute.
- Press the "MEMORY" button to read the next line of stored data. The data called up on the display is numbered accordingly.
- To discontinue reading the data, press "START/STOP" button to normal time.

---

## Maintenance.

- Use a piece of cloth with water or mild cleansing agent to clean the case and then use a piece of dry cloth to wipe it dry. Use a piece of dry cloth to wipe the cuff when it is dirty
- Do not use gas or any strong cleansers in cleaning.
- When the unit is not to be used for a long time, remove the batteries. (Leaking of battery liquid can cause trouble.)

## Warning

- Do not in any way twist the arm cuff.
- Do not start the start measurement when the arm cuff is not wrapped around the arm.
- Do not attempt to disassemble or change any parts of the monitor including the cuff.
- Do not drop the product or put it through strong impact.
- The cuff size is suitable for 23-33 cm / 9-13 inches adult only.
- The device is not supposed to be used when your arm has any wound or injury.
- In case the cuff kept pumping up non-stop, open the cuff at once.
- We do not recommend user to unpack the device by himself/herself to substitute of a component different from that supplied might result in measurement error.
- If service is needed please contact your local Seagull reseller.
- This product is not suitable for pregnant women.
- To avoid accidental strangulation, keep this product away from children and do not drape tube around neck.

---

**Recalibration Notice :**

To ensure precise measurement reading, recalibration of the device is recommended after 2 years from the date of manufacturing. Shipping plus handling cost and recalibration service fee shall be charged accordingly. Please contact us at (+45) 5621 6608.

---

Method of measurement.....	Oscillometric
Model Number.....	HL888HA
Range of measurement.....	Pressure 0-300 mmHg
Pulse .....	40 - 199 beats/minute
Accuracy.....	Pressure +/- 3 mmHg
Pulse .....	+/- 5% Max
Inflation.....	Deluxe Automatic
Exhaust.....	Automatic exhaust valve
Display.....	Systolic, Diastolic, Pulse rate
Sets of memory.....	99 sets
Cuff size.....	Arm circumference approx.
.....	23-33 cm/9-13 inch
Operating temperature.....	+10°C - +40°C, less than 85% R.H.
Storage temperature.....	-20°C - +70°C, less than 85% R.H.
Power Supply.....	6 VDC, 4xAA (1,5V) Alkaline
.....	batteries
Battery Life.....	300 measurements
Accessory.....	4 batteries, gift box, manual,
.....	storage bag

**BF**

**Classification:**

- Internally powered equipment
- BF type applied part
- IPX0
- Not suitable for use in presence of flammable anesthetic mixture with air or with Oxygen or nitrous oxide
- Continuous operation with short-time loading

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipments, do not use the device near a mobile phone or microwave oven.

Discard the used product to the recycling collection point according to local regulations.

**Manufacturer:**

HEALTH & LIFE Co., Ltd.

9F, No. 186, Jian Yi Road, Zhonghe District, 23553  
New Taipei City, Taiwan

---

**EC REP**

**EMERGO EUROPE**

Molenstraat 15, 2513 BH, The Hague, The  
Netherlands

TEL : +31-70-3458570

FAX : +31-70-3467299

BF

**Classification:**

- Internally powered equipment
- BF type applied part
- IPX0
- Not suitable for use in presence of flammable anesthetic mixture with air or with Oxygen or nitrous oxide
- Continuous operation with short-time loading

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipments, do not use the device near a mobile phone or microwave oven.

Discard the used product to the recycling collection point according to local regulations.

**Manufacturer:**

HEALTH & LIFE Co., Ltd.

9F, No. 186, Jian Yi Road, Zhonghe District, 23553  
New Taipei City, Taiwan

---

# Appendix

<b>Guidance and manufacturer's declaration - electromagnetic emissions</b>			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>	
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Not applicable		
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable		
<b>Guidance and manufacturer's declaration - electromagnetic immunity</b>			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact  ±8 kV air	±6 kV contact  ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

### Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. <b>Recommended separation distance</b> $d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz  $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a)</sup> should be less than the compliance level in each frequency range. <sup>b)</sup> Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a). Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b). Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### **Recommended separation distances between portable and mobile RF communications equipment and the device**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

---

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# Blood Pressure Diary

---

Date	Time	Systolic/Diastolic	Pulse
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

# Blood Pressure Diary

---

Date	Time	Systolic/Diastolic	Pulse
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

# Blood Pressure Diary

---

Date	Time	Systolic/Diastolic	Pulse
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	
		/	

Distributor: SEAGULL Aps  
Theilgaards Allé 3c  
DK-4600 Køge  
Denmark  
+45 5621 6608

[www.seagull-healthcare.dk](http://www.seagull-healthcare.dk)