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Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100

Directions for use

901042 Automated Blood Pressure System, Software Version 1.0

© 2016 Welch Allyn. All rights are reserved. To support the intended use of the product described in this publication, the purchaser of the product is permitted to copy this publication, for internal distribution only, from the media provided by Welch Allyn. No other use, reproduction, or distribution of this publication, or any part of it, is permitted without written permission from Welch Allyn.

Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

The *Bluetooth*® word mark and logos are registered trademarks owned by *Bluetooth* SIG, Inc. and any use of such marks by Welch Allyn is under license.

Software in this product is Copyright 2016 Welch Allyn or its vendors. All rights are reserved. The software is protected by United States of America copyright laws and international treaty provisions applicable worldwide. Under such laws, the licensee is entitled to use the copy of the software incorporated with this instrument as intended in the operation of the product in which it is embedded. The software may not be copied, decompiled, reverse-engineered, disassembled, or otherwise reduced to human-perceivable form. This is not a sale of the software or any copy of the software; all right, title, and ownership of the software remain with Welch Allyn or its vendors.

For information about any Welch Allyn product, contact your local Welch Allyn representative: www.welchallyn.com/about/company/locations.htm.

723807

This manual applies to **REF** 901042 Automated Blood Pressure System.

DIR 80020549 Ver. B Revised 2016-01
 Manufactured for Welch Allyn, Inc.
 4341 State Street Road
 Skaneateles Falls, NY 13153 USA

Guangdong Transtek Medical Electronics Co., Ltd.
 Zone A, No. 105, Dongli Road
 Torch Development District
 Zhongshan, 528437, Guangdong, China

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Introductory information

Introduction

The Welch Allyn remote monitoring blood pressure device is indicated for adult use in the home or domestic setting only.

The Welch Allyn blood pressure monitor features blood pressure measurement, pulse rate measurement, and data management using *Bluetooth*® Smart technology.

The design provides you with two years of reliable service. Readings taken by the monitor are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method.

This *Directions for use* contains important safety and care information and provides step by step instructions for using the monitor. Read the manual thoroughly before using the monitor.

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Indications for use

The Welch Allyn remote monitoring blood pressure device is a digital monitor intended for use in measuring blood pressure and heart rate in adult patient population with arm circumferences between 22 cm to 42 cm (approximately 9 to 17 inches).

The device detects the appearance of irregular heartbeats during measurement and gives a warning signal with readings.

Contraindications

This device is contraindicated for any person who is connected to a wearable or implantable electronic device or instrument, such as a pacemaker or defibrillator.

This blood pressure monitor is not intended to be a diagnostic device. Contact your physician if hypertensive values are indicated.

Symbols

Documentation symbols

Warning: The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death.

Caution: The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.

Follow instructions/directions for use (DFU) -- mandatory action.

Power symbols

Direct current

Shipping, storing, and environment symbols

Relative humidity limits

Recycle the product separate from other disposables

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Cuff symbols

 ARTERY	Artery symbol
	Range
INDEX 	Index edge
	Limb circumference (Minimum/Maximum)
LOT	Lot code
	Not made with natural rubber latex

Miscellaneous symbols

	Type BF applied parts		Serial number
REF	Product Identifier		Manufacturer
#	Reorder number, Model number		Non-ionizing electromagnetic radiation
GTIN	Global Trade Identification Number	IP22	The ingress protection: the device could protected against solid foreign objects of 12.5mm and greater, and against vertically falling water drops when ENCLOSURE tilted up to 15°
	Class II equipment		

About warnings and cautions

Caution statements can appear on the Welch Allyn blood pressure device, the packaging, the shipping container, or in this *Directions for use*.

The Welch Allyn blood pressure device is safe for patients and clinicians when used in accordance with the instructions and caution statements presented in this *Directions for use*.

Before using the device, you must familiarize yourself with all cautions, with the steps to power up the device, and with the sections of this *Directions for use* that pertain to your use of the device. In addition to reviewing the general cautions presented in the next section, you must also review the more specific cautions that appear throughout the manual in conjunction with setup/startup, operation, and maintenance tasks.

- Failure to understand and observe any warning statement in this manual could lead to patient injury or illness.
- Failure to understand and observe any caution statement in this manual could lead to damage to the equipment or other property, or loss of patient data.

General warnings and cautions

Warning Patient injury risk. The device is not suitable for measuring the blood pressure of neonatal infants or children. Ask your doctor before using it with older children.

Warning Patient injury risk. The device is not suitable for use on pregnant women, patients with implanted electrical devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral arterial disease, and patients undergoing intravascular therapy or arterio-venous shunt or people who received a mastectomy. Consult your doctor prior to using the unit if you suffer from any of these illnesses or conditions.

Warning Patient injury risk. Do not take any therapeutic measures on the basis of a self measurement. Never alter the dose of a medicine prescribed by a doctor. Consult your doctor if you have any question about your blood pressure.

Warning Patient injury risk. Keep the unit out of reach of infants, children, or pets, since inhalation or swallowing of small parts is dangerous or even fatal.

Warning Patient injury risk. The device is not intended to be a self-measurement station used in a public place or center.

Warning Patient injury risk. Do not burn batteries. Batteries may leak or explode.

Warning Patient injury risk. If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the Standby symbol immediately to release the air from the cuff. Loosen and remove the cuff from your arm.

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-  **Warning** Patient injury risk. On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure applied to the arm (cuff pressure >300mmHg or constant pressure >15mmHg for more than 3 minutes) may lead to an ecchymosis.
-  **Warning** Patient injury risk. This unit is not suitable for continuous monitoring during medical emergencies or operations
-  **Warning** Patient injury risk. Too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.
-  **Warning** Patient injury risk. Do not kink the connection tube during use. The cuff pressure may continuously increase which can prevent blood flow and result in harmful injury.
-  **Warning** Patient injury risk. Never apply the cuff over damaged skin. This can cause further injury.
-  **Warning** Patient injury risk. Do not touch output of the batteries/adapter and the patient simultaneously.
-  **Warning** Patient injury risk. Do not use the unit in case of existing polyester resp. synthetic allergies.
-  **Warning** Patient injury risk. Excessive cable and hose lengths could cause strangulation due to excessive length.
-  **Warning** Patient injury risk. Do not connect the air hose to other medical equipment. This could cause air to be pumped into intravascular systems or high pressure, what could lead to dangerous injuries.
-  **Warning** Patient injury risk. Dispose of accessories, detachable parts, and the ME equipment according to the local guidelines.
-  **Warning** Patient injury risk. Do not service or perform any maintenance while the device is in use.
-  **Warning** Patient injury risk. Use accessories and detachable parts specified/ authorized by the manufacturer. Otherwise, it may cause damage to the unit or danger to the user/patients.
-  **CAUTION** This device is intended only for adult use in homes.
-  **CAUTION** This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than arm or for functions other than obtaining a blood pressure measurement.
-  **CAUTION** Do not wrap the cuff on the same arm which other monitoring ME equipment is applied simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring ME equipment..
-  **CAUTION** The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air of oxygen or nitrous oxide.

-  **CAUTION** This device cannot be used with HF surgical equipment at the same time.
-  **CAUTION** To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
-  **CAUTION** Manufacturer will make available on request circuit diagrams, component parts list etc.
-  **CAUTION** The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.
-  **CAUTION** Before use, make sure the device functions safely and is in proper working condition.
-  **CAUTION** Use the device in the environment which was provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.
-  **CAUTION** Do not attempt to repair the unit yourself in the event of malfunctions. Only have repairs carried out by authorized service centers.
-  **CAUTION** Report any unexpected operation or events to the manufacturer.
-  **CAUTION** The device does not need to be calibrated in two years of reliable service.
-  **CAUTION** Use a soft cloth to clean the entire unit. Do not use any abrasive or volatile cleaners.

Content list

The following items are in the device box:

- Blood pressure monitor
- Cuff
- (4) AA alkaline batteries
- *Directions for use*

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

Controls and indicators

Monitor front

No.	Feature	Description
1	Cuff (Type BF applied part)	Attaches to arm of patient
2	Air hose	Connects the cuff to the air connector plug
3	Air connector plug	Air hose connects to the air connector plug
4	Power button	Powers on the monitor, starts the blood pressure measurement, and – when held – initiates a <i>Bluetooth</i> pairing.
5	LCD Display	Displays blood pressure reading and other pertinent information regarding the reading

Monitor back

No.	Feature	Description
1	Direct current power connection	Connects the monitor to a mains power outlet (not included)
2	Battery compartment (behind cover)	Houses 4 AA alkaline batteries

Power options

CAUTION In order to get the best effect and protect your monitor, please use the correct battery or Welch Allyn specified power adaptor.

The monitor is powered by one of two sources:

- 4 AA alkaline batteries
- AC adaptor powered mode (6v 1A) (not included)

Insert or replace the batteries

CAUTION Remove the batteries if the device is not to be used for some time.

CAUTION Dispose of old batteries by following your local recycling guidelines.

CAUTION Do not burn batteries. Batteries may leak or explode.

If you are not using the accessory power cord, you must install 4 AA alkaline batteries before the device can be used.

Replace the batteries when any of the following occurs.

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
- The display does not light up
- Slide off the battery cover.
 - Install the batteries by matching the correct polarity as shown in the diagram.

- Replace the cover.

Screen elements

The liquid crystal display (LCD) may indicate any of the following: systolic blood pressure (mmHg), diastolic blood pressure (mmHg), pulse rate (bpm), and battery charge level.

Symbol	Description	Explanation
	Systolic blood pressure	High pressure result mmHG = Measurement Unit of the blood pressure
	Diastolic blood pressure	Low pressure result mmHG = Measurement Unit of the blood pressure

Symbol	Description	Explanation
	Pulse	Pulse/Min
	Irregular heartbeat	Irregular heartbeat
	Heart beat	Monitor is detecting a heartbeat during measurement
	Battery indicator	Indicates the current battery charge
	Low battery indicator	Indicates the current battery charge
	Data transmission error	Data transmission error
	Motion indicator	Motion may result in an inaccurate measurement.
	Out of range	Either SYS > 230mmHg or DIA > 130mmHg. The symbol may appear in either the SYS or DIA area of the screen. Alarm priority = Low

Segments	Pression
1	>=0
2	>=40
3	>=80
4	>=120
5	>=140
6	>=160
7	>=180
8	>=200

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Use the monitor

Startup

The patient is an intended operator of the device. The patient can measure, transmit data, change the battery under normal circumstances, and maintain the device and its accessories according to the user manual.

Press the **Power** button to turn the device on or off. Upon each power up, the display lights up and displays "888."

If a system error is detected, the device displays a system error message beginning with the letter "E" and followed by a system fault code to assist service and engineers in diagnosing the problem.

Pair a communication hub with the monitor

With advanced Bluetooth 4.0 Technology applied, a communication hub is capable of receiving your personal health data. Contact your Remote Patient Monitoring Solution provider for connectivity instructions and troubleshooting.

1. Confirm that both the mobile or portable equipment and the monitor are powered on.

Note *Bluetooth* wireless connectivity is only available for devices which support *Bluetooth* 4.0 technology.

2. Press and hold the **Power/Standby** button for 2 seconds.

The *Bluetooth* pairing sequence is activated. The dotted line surrounding the *Bluetooth* symbol flashes while discovering the monitor.

If successful, a steady  symbol surrounded by a dotted line appears on the LCD screen. If unsuccessful, the  appears.

Blood pressure items

Identify and have available the blood pressure monitor, the blood pressure cuff, and the blood pressure air hose.

Attach the blood pressure cuff

The preferred blood pressure measurement site is the upper arm. Keep your arm relaxed and motion-free during measurement(s).

Complete these steps to attach the blood pressure cuff to your arm.

1. Bare your arm to attach the cuff.
2. Sit comfortably with your legs uncrossed, feet flat on the floor, and back and arm supported.

The center of the cuff should be at the same level as the right atrium of your heart.

3. Attach the cuff to your upper arm with the air hose off-center toward the inner side of your arm in line with your little finger.

The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.

4. Sit comfortably with your arm resting on a flat surface.
5. Rest 5 minutes before beginning the measurement.
6. Wait at least 3 minutes between measurements. This allows your blood circulation to recover.

For a meaningful comparison, try to take measurements under consistent conditions. For example, take daily measurements at approximately the same time, on the same upper arm, or as directed by your physician.

Start the measurement

The preferred blood pressure measurement site is the upper arm. Keep the patient's arm relaxed and motion-free during measurement(s).

When the monitor is off, press the **Power/Standby** button to power up the monitor and begin the measurement.

It will finish the entire measurement. During inflation, the heart icon in the upper-left corner blinks in keeping with your pulse rate.

After successful transfer, the device powers off the *Bluetooth* radio and the icon (and rectangle) are removed.

If you press and release the **Power/Standby** button, the device powers down.

If there is no operation, after a 10 seconds of inactivity, the device powers down.

Additionally, during inflation the progress meter to the right of the digits builds vertically as the pressure increases according to the table below.

Segments	Pressure
1	>=0
2	>=40
3	>=80
4	>=120

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If you press the **Power/Standby** button to stop the measurement during the measurement, the SYS and DIA numbers are cleared and replaced with dotted lines.

CAUTION The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other record numbers are increased by one digit (for example, 2 becomes 3, and so on), and the oldest record (99) is dropped from the list.

Blood pressure notes

Activities that could cause inaccurate measurements

- Taking a measurement within 1 hour of eating or drinking
- Taking a measurement within 20 minutes of taking a bath
- Being in a very cold environment
- Taking a measurement immediately after drinking tea, coffee, or smoking
- Taking a measurement while talking or moving your fingers
- Taking a measurement when you need to urinate

What you need to pay attention to when you measure your blood pressure at home

- If the cuff is attached properly.
- If the cuff is too tight or too loose.
- If the cuff is attached on the upper arm.
- If you feel anxious, take 2–3 deep breaths before beginning. Relax for 4–5 minutes until you calm down. This will result in a more accurate measurement.

What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of your heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

What is the standard blood pressure classification?

The chart below is the standard blood pressure classification published by American Heart Association (AHA).

Blood Pressure Category	Systolic mmHg (upper#)		Diastolic mmHg (lower#)
Normal	less than 120	and	less than 80
Prehypertension	120–139	or	80–89
High blood pressure (Hypertension) Stage 1	140–159	or	90–99
High blood pressure (Hypertension) Stage 2	160 or higher	or	100 or higher
Hypertensive Crisis (Emergency care needed)	Higher than 180	or	Higher than 110

CAUTION Only a physician can tell your normal BP range. Contact a physician if your measuring result falls out of the range. Only a physician can tell whether your blood pressure value has reached a dangerous point.

Irregular Heartbeat Detector

An irregular heartbeat (IHB) is detected when a heartbeat rhythm varies while the unit is measuring the systolic and diastolic blood pressure. During each measurement, the monitor records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 15, the irregular heartbeat symbol appears on the symbol when the measurement results are displayed.

CAUTION The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. The device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies multiple times everyday. It is also affected by the way you tie your cuff and your measurement position, so please take the measurement under the same conditions.
2. If the person takes medicine, the pressure will vary more.
3. Wait at least 3 minutes for another measurement.

Why do I get a different blood pressure at home compared to the hospital?

The blood pressure is different even throughout the day due to weather, emotion, exercise, and the like. Also, there is the "white coat" effect, which means blood pressure may be increased in clinical settings.

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Is the result the same if measuring on the right arm?

It is alright for both arms, but there will be some different results for different people. We suggest you measure the same arm every time.

Maintenance and service

Maximize performance

CAUTION Use a soft cloth to clean the entire unit. Do not use any abrasive or volatile cleaners.

In order to get the best performance, please follow the steps below.

- Put in a dry place and avoid the sun.
- Avoid touching water. Clean the monitor with a dry cloth.
- Avoid intense shaking and collisions.
- Avoid dusty and unstable temperature environments.
- Do not use wet cloths to remove dirt.
- Avoid washing the cuff.

Clean the monitor

Warning Patient injury risk. Keep the device clean. Clean the device after each use. Contact with contaminated equipment could spread infection.

Routinely clean the monitor according to the facility's protocols and standards or local regulations.

The following agents are compatible with the device:

- CaviWipes or Sani-Cloth Plus
- 70% isopropyl alcohol
- 10% chlorine bleach/90% water solution

Troubleshooting

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the device is not operating as you think it should, check here before arranging for servicing.

Problems and error messages

Problem	Symptom	Root cause	Solution
No power	Display will not light up	Batteries are exhausted.	Replace with new batteries
		Batteries are inserted incorrectly.	Insert the batteries correctly


Problem	Symptom	Root cause	Solution
		AC adaptor is inserted incorrectly.	Insert the AC adaptor tightly
Low batteries	The display indicates the "0-LO" message, pauses for 3 seconds. The battery icon shows empty (does not flash.)	Batteries are low.	Replace with new batteries
Error messages	E 01 shows	The cuff is not secure	Readjust the cuff, relax for a moment, and then measure again
	E 02 shows	The cuff is very tight	Refasten the cuff and then measure again
	E 03 shows	There is too much pressure in the cuff	Refasten the cuff and then measure again
	E 10 or E 11 shows	The monitor detected motion while measuring	Relax for a moment and then measure again
	E 20 shows	The measurement process does not detect a pulse signal	Loosen the clothing on the arm and then measure again
	E 21 shows	Measurement incorrect	Relax for a moment and then measure again

Alarm indications

Symptom	Root cause	Solution
	Out of range. Either SYS > 230mmHg or DIA > 130mmHg. The symbol may appear in either the SYS or DIA area of the screen.	Press and hold the Power button to reestablish the Bluetooth connectivity. Measure again. If the problem persists, contact your physician. Alarm priority = Low !

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Symptom	Root cause	Solution
	Out of range. Either SYS < 60mmHg or DIA < 40mmHg. The symbol may appear in either the SYS or DIA area of the screen.	Press and hold the Power button to reestablish the Bluetooth connectivity. Measure again. If the problem persists, contact your physician. Alarm priority = Low !

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Appendices

Specifications

Warning: No modification of this equipment is allowed.

Item	Specification
Power supply: Battery powered mode:	6VDC 4*AA batteries
Power supply: AC adaptor powered mode:	Input: 100-240V, 50-60Hz, 400mA; Output: 6V, 1A
Display mode	Digital LCD V.A 68mm*90mm
Measurement mode	Oscillometric testing mode
Measurement range	Rated cuff pressure: 0mmHg-300mmHg (0kPa-40kPa) Measurement pressure: SYS: 60-230mmHg DIA: 40-130mmHg Pulse value: (40-199)beat/minute
Accuracy	Pressure: 5°C-40°C within±0.4kPa(3mmHg) Pulse value:±4%
Normal working condition	Temperature:5°C to 40°C Relative humidity ≤ 85%RH Atmospheric pressure: 86kPa to 106kPa
Storage & transportation condition	Temperature:-20°Cto 60°C Relative Humidity: 10% RH-93% RH Atmospheric Pressure: 50kPa to 106 kPa
Measurement perimeter of the upper arm	About 22cm-42cm

Item	Specification
Net Weight	Approx.300g (Excluding the dry cells)
External dimensions	Approx.92mm*140mm*46mm
Attachment	4*AA alkaline batteries, <i>Directions for use</i>
Degree of protection	Type BF applied part
Protection against ingress of water	IP22
Software Version	Version 1.0

General radio compliance

Item	Specification
Bluetooth Module No.	AW8001
RF Frequency Range	2402 MHz to 2480 MHz
Output Power Range	-2.27dBm to -5.88 dBm
Supply Voltage	1.8-3.6 V
Antenna Gain	1.0dBi
Transmitting Distance	10 meters, (30 feet)

Federal Communication Commission (FCC) Interference Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions.

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

This equipment (FCC ID: OU9TMB1591-B) has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures.

- 1.Reorient or relocate the receiving antenna.
- 2.Increase the separation between the equipment and receiver.

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Bluetooth activé

3. Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

FCC Radiation Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Industry Canada (IC) compliance

To ensure compliance with FCC and Industry Canada RF exposure requirements, this device must be installed in a location where the antennas of the device will have a minimum distance of at least 20 cm from all persons. Using higher gain antennas and types of antennas not certified for use with this product is not allowed. The device shall not be co-located with another transmitter.

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This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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This radio transmitter (IC ID: 12725A-TMB1591B) has been approved by Industry Canada to operate with the antenna listed in the specification table.

This Class B digital apparatus complies with Canadian ICES-003.

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Complied standards list

Item	Standard
Risk management	ISO/EN 14971 Medical devices — Application of risk management to medical devices
Labeling	ISO/EN 15223-1 Medical devices. Symbols to be used with medical device labels, labeling and information to be supplied. General requirements
User manual	EN 1041 Medical equipment manufacturers to provide information

Item	Standard
General Requirements for Safety	IEC 60601-1+A1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance IEC 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic compatibility	IEC/EN 60601-1-2 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
Performance requirements and clinical investigation	IEC 80601-2-30 Medical electrical equipment- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers ISO81060-2 Non-invasive sphygmomanometers — Part 2: Clinical validation of automated measurement type
Software life-cycle processes	IEC/EN 62304+AC: Medical device software - Software life cycle processes
Usability	IEC 62366 Medical devices - Application of usability engineering to medical devices (IEC 62366) IEC 60601-1-6 Medical electrical equipment - Part 1 -6: General requirements for basic safety and essential performance - collateral standard : Usability

Emissions and immunity information

Electromagnetic emissions

The Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 is intended for use in the electromagnetic environment specified below. The customer or user of the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 should assure that it is used in such an environment.

This equipment needs to be installed and put into service in accordance with the information provided in the user manual.

CAUTION Wireless communications equipment, such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d=3,3m$ away from the equipment.

Note As indicated in Table 6 of IEC 60601-1-2:2007 for ME Equipment, a typical cell phone with a maximum output power of 2 W yields $d=3,3m$ at an IMMUNITY LEVEL of 3V/m).

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 uses RF energy only for its internal function. Therefore, its RF emissions

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Electromagnetic emissions

are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11	Class B	The Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 is suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Electromagnetic immunity

The Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 is intended for use in the electromagnetic environment specified below. The customer or the user of the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	± 1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 seconds	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 requires continued operation during power mains interruptions, it is recommended that the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 be powered from an uninterruptible power supply or a battery.

Electromagnetic immunity

Power frequency (50Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
---	------	------	---

Note: U_1 is the AC mains voltage prior to application of the test level.

Electromagnetic immunity

The Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 is intended for use in the electromagnetic environment specified below. The customer or the user of the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
---------------	----------------------	------------------	--

Portable and mobile RF communications equipment should be used no closer to any part of the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Recommended separation distance

Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = (1.167) \sqrt{P}$
----------------------------	-----------------------------	--------	------------------------

Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 1 GHz	3 V/m	$d = (1.167) \sqrt{P}$ 80 to 800 MHz
---------------------------	-----------------------	-------	--------------------------------------

$d = (2.333) \sqrt{P}$ 800 MHz to 2,5 GHz

where P is the maximum output power rating of the transmitter in watts (W) and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 is used exceeds the applicable RF compliance level above, the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 should be observed to verify normal operation. If abnormal performance is observed,

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4341 State Street Road

Skaneateles Falls, NY 13153 États-Unis

Guangdong Transtek Medical Electronics Co., Ltd.
Zone A, No. 105, Dongli Road
Torch Development District
Zhongshan, 528437, Guangdong, Chine

www.welchallyn.com

WelchAllyn
Advancing Frontline Care™

Electromagnetic immunity

additional measures may be necessary, such as reorienting or relocating the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100.

^bOver the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100

The Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Welch Allyn Remote Monitoring Blood Pressure Device RPM-BP100 as recommended below, according to the maximum output power of the communications equipment.

Separation distance according to frequency of transmitter (m)

Rated max. output power of transmitter (W)	150 kHz to 80 MHz $d = (1.167) \sqrt{P}$	80 MHz to 800 MHz $d = (1.167) \sqrt{P}$	800 MHz to 2.5 GHz $d = (2.333) \sqrt{P}$
0.01	0.167	0.167	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.690	3.690	7.338
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

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Warranty

Welch Allyn will warranty the Transtek blood pressure monitor to be free of defects in material and workmanship and to perform in accordance with manufacturer specifications for the period of two years from the date of purchase from Welch Allyn or its authorized distributors or agents.

The warranty period shall start on the date of purchase. The date of purchase is: 1) the invoiced ship date if the device was purchased directly from Welch Allyn, 2) the date specified during product registration, 3) the date of purchase of the product from a Welch Allyn authorized distributor as documented from a receipt from said distributor.

This warranty does not cover damage caused by: 1) handling during shipping, 2) use or maintenance contrary to labeled instructions, 3) alteration or repair by anyone not authorized by Welch Allyn, and 4) accidents.

The product warranty is also subject to the following terms and limitations.

- Accessories are not covered by the warranty. Refer to the *Directions for use* provided with individual accessories for warranty information.
- Shipping cost to return a device to a Welch Allyn service center is not included.
- A service notification number must be obtained from Welch Allyn prior to returning any products or accessories to Welch Allyn's designated service centers for repair. To obtain a service notification number, contact Welch Allyn Technical Support at www.welchallyn.com/about/company/locations.htm.

Approved accessories

Item	Description
RPM-BPACC-02	D-Ring Standard Wide Cuff (22-42cm)
RPM-BPACC-04	RPM BP AC Adapter. This adapter is an alternate power source for the blood pressure monitoring device. (Model number UE08WCP-06100SPA)

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